



***In vitro* Inflammation Model induced by Cytokines**

Reference:

Enrriquez-de-Salamanca A., Calder V., Gao J., Galatowicz G., García-Vázquez C., Fernández I., Stern M.E., Diebold Y., Calonge M. Cytokine responses by conjunctival epithelial cells: An in vitro model of ocular inflammation. *Cytokine*, 2008; 44: 160-167.

Fundamentals:

The ***In vitro* Inflammation Model induced by Cytokines** allows identifying the presence of an inflammatory process in cell cultures through an increased expression of proinflammatory metabolites such as IL-6, following exposure to cytokines such as TNF- α .

Technique:

Seed Human Corneal Epithelial (HCE) or Normal Human Conjunctiva (IOBA-NHC) cells in a 24-well plate (80,000 cells/well in 500 μ l culture medium/well). Gently shake plates after seeding in order to homogenize cell distribution over the well surface. Let them grow for 48 h at standard conditions (37 $^{\circ}$ C, 5% CO₂/95% air atmosphere).

Replace HCE or IOBA-NHC culture media for serum-free, nonsupplemented medium (DMEM/F-12). Maintain cells in the incubator for 24 h.

Afterwards, set experimental conditions:

Positive control: TNF- α in serum-free, nonsupplemented medium (25 ng/ml). Use a handy TNF- α concentration in mother solution like 100 μ g/ml to prepare the working solution.

Negative control: serum-free, nonsupplemented medium.


Samples: in serum-free, nonsupplemented medium.

Maintain cells in the incubator for 24 h.

Collect cell culture supernatants and maintain at -80° C until use.

Quantify IL-6 secretion using the human IL-6 ELISA Kit.

Timeline:



3D-NET (EU FP7 2013 - 612218/3D-NET)

3DNET SOPs

HOW-TO-DO PRACTICAL GUIDE



Note: HCE and IOBA-NHC media are described in “HowToDo-human ocular cell lines”.